

DODA Antibody



PACO58220

Product Information

Size:

50ug

Reactivity:

Portulaca grandiflora

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Opens the cyclic ring of dihydroxy-phenylalanine (DOPA) between carbons 4 and 5, thus producing an unstable seco-DOPA that rearranges nonenzymatically to betalamic acid, .

Gene ID:

DODA

Uniprot

Q7XA48

Synonyms:

4,5-DOPA dioxygenase extradiol (EC 1.13.11.29), DODA

Immunogen:

Recombinant Portulaca grandiflora 4,5-DOPA dioxygenase extradiol protein (1-271AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A