

SETD3 Antibody, HRP conjugated



PACO58229

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Histone methyltransferase that methylates 'Lys-4' and 'Lys-36' of histone H3 (H3K4me and H3K36me). Acts as a transcriptional activator. Plays an important role in the transcriptional regulation of muscle cell differentiation via interaction with MYOD1.

Gene ID:

SETD3

Uniprot

Q86TU7

Synonyms:

Histone-lysine N-methyltransferase setd3 (EC 2.1.1.43) (SET domain-containing protein 3), SETD3, C14orf154

Immunogen:

Recombinant Human Histone-lysine N-methyltransferase setd3 protein (166-265AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A