

KTN1 Antibody, FITC conjugated



PACO58238

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for kinesin thus involved in kinesin-driven vesicle motility. Accumulates in integrin-based adhesion complexes (IAC) upon integrin aggregation by fibronectin.

Gene ID:

KTN1

Uniprot

Q86UP2

Synonyms:

Kinectin (CG-1 antigen) (Kinesin receptor), KTN1, CG1 KIAA0004

Immunogen:

Recombinant Human Kinectin protein (978-1116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A