

# PNPLA5 Antibody, FITC conjugated



PACO58250

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Lipid hydrolase.

**Gene ID:**

PNPLA5

**Uniprot**

Q7Z6Z6

**Synonyms:**

Patatin-like phospholipase domain-containing protein 5 (EC 3.1.1) (GS2-like protein), PNPLA5, GS2L

**Immunogen:**

Recombinant Human Patatin-like phospholipase domain-containing protein 5 protein (90-205AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A