

PNPLA5 Antibody, Biotin conjugated



PACO58251

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Lipid hydrolase.

Gene ID:

PNPLA5

Uniprot

Q7Z6Z6

Synonyms:

Patatin-like phospholipase domain-containing protein 5 (EC 3.1.1) (GS2-like protein), PNPLA5, GS2L

Immunogen:

Recombinant Human Patatin-like phospholipase domain-containing protein 5 protein (90-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A