

CSMD3 Antibody, FITC conjugated



PACO58258

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in dendrite development.

Gene ID:

CSMD3

Uniprot

Q7Z407

Synonyms:

CUB and sushi domain-containing protein 3 (CUB and sushi multiple domains protein 3), CSMD3, KIAA1894

Immunogen:

Recombinant Human CUB and sushi domain-containing protein 3 protein (233-380AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A