

# SAT2 Antibody, Biotin conjugated



PACO58311

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Enzyme which catalyzes the acetylation of polyamines. Substrate specificity: norspermidine > spermidine = spermine >> N(1)acetylspermine = putrescine.

**Gene ID:**

SAT2

**Uniprot**

Q96F10

**Synonyms:**

Diamine acetyltransferase 2 (EC 2.3.1.57) (Polyamine N-acetyltransferase 2) (Spermidine/spermine N(1)-acetyltransferase 2) (Thialysine N-epsilon-acetyltransferase), SAT2, SSAT2

**Immunogen:**

Recombinant Human Diamine acetyltransferase 2 protein (1-170AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A