

SLC5A8 Antibody, FITC conjugated



PAC058322

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SLC5A8

Uniprot

Q8N695

Synonyms:

Sodium-coupled monocarboxylate transporter 1 (Apical iodide transporter) (Electrogenic sodium monocarboxylate cotransporter) (Sodium iodide-related cotransporter) (Solute carrier family 5 member 8), SLC5A8, AIT SMCT SMCT1

Immunogen:

Recombinant Human Sodium-coupled monocarboxylate transporter 1 protein (540-610AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A