

PANX2 Antibody, Biotin conjugated



PACO58339

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Structural component of the gap junctions and the hemichannels.

Gene ID:

PANX2

Uniprot

Q96RD6

Synonyms:

Pannexin-2, PANX2

Immunogen:

Recombinant Human Pannexin-2 protein (370-470AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A