

KCNV2 Antibody, FITC conjugated



PACO58386

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potassium channel subunit. Modulates channel activity by shifting the threshold and the half-maximal activation to more negative values.

Gene ID:

KCNV2

Uniprot

Q8TDN2

Synonyms:

Potassium voltage-gated channel subfamily V member 2 (Voltage-gated potassium channel subunit Kv8.2), KCNV2

Immunogen:

Recombinant Human Potassium voltage-gated channel subfamily V member 2 protein (1-155AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A