

LIPH Antibody, HRP conjugated



PACO58393

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Hydrolyzes specifically phosphatidic acid, (PA) to produce 2-acyl lysophosphatidic acid, (LPA; a potent bioactive lipid mediator) and fatty acid, Does not hydrolyze other phospholipids, like phosphatidylserine (PS), phosphatidylcholine (PC) and phosphatidylethanolamine (PE) or triacylglycerol (TG).

Gene ID:

LIPH

Uniprot

Q8WWY8

Synonyms:

Lipase member H (LIPH) (EC 3.1.1) (LPD lipase-related protein) (Membrane-associated phosphatidic acid, selective phospholipase A1-alpha) (mPA-PLA1 alpha) (Phospholipase A1 member B), LIPH, LPDLR MPAPLA1 PLA1B

Immunogen:

Recombinant Human Lipase member H protein (175-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A