

CLPTM1L Antibody, FITC conjugated



PACO58410

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Enhances cisplatin-mediated apoptosis, when overexpressed.

Gene ID:

CLPTM1L

Uniprot

Q96KA5

Synonyms:

Cleft lip and palate transmembrane protein 1-like protein (CLPTM1-like protein) (Cisplatin resistance-related protein 9) (CRR9p), CLPTM1L, CRR9

Immunogen:

Recombinant Human Cleft lip and palate transmembrane protein 1-like protein (32-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A