

# NKD1 Antibody, Biotin conjugated



PACO58419

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Cell autonomous antagonist of the canonical Wnt signaling pathway. May activate a second Wnt signaling pathway that controls planar cell polarity.

**Gene ID:**

NKD1

**Uniprot**

Q969G9

**Synonyms:**

Protein naked cuticle homolog 1 (Naked-1) (hNkd) (hNkd1), NKD1, NKD

**Immunogen:**

Recombinant Human Protein naked cuticle homolog 1 protein (269-384AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A