

SLC34A3 Antibody, HRP conjugated



PACO58425

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in actively transporting phosphate into cells via Na(+) cotransport in the renal brush border membrane. Probably mediates 20-30% of the apical influx.

Gene ID:

SLC34A3

Uniprot

Q8N130

Synonyms:

Sodium-dependent phosphate transport protein 2C (Sodium-phosphate transport protein 2C) (Na(+)-dependent phosphate cotransporter 2C) (Sodium/inorganic phosphate cotransporter IIC) (Sodium/phosphate cotransporter 2C) (Na+)/Pi cotransporter 2C) (NaPi-2c) (Solute carrier family 34 member 3), SLC34A3, NPT2C NPTIIC

Immunogen:

Recombinant Human Sodium-dependent phosphate transport protein 2C protein (1-76AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A