

# SLC34A3 Antibody, Biotin conjugated



PACO58427

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in actively transporting phosphate into cells via Na(+) cotransport in the renal brush border membrane. Probably mediates 20-30% of the apical influx.

**Gene ID:**

SLC34A3

**Uniprot**

Q8N130

**Synonyms:**

Sodium-dependent phosphate transport protein 2C (Sodium-phosphate transport protein 2C) (Na(+)-dependent phosphate cotransporter 2C) (Sodium/inorganic phosphate cotransporter IIC) (Sodium/phosphate cotransporter 2C) (Na+)/Pi cotransporter 2C) (NaPi-2c) (Solute carrier family 34 member 3), SLC34A3, NPT2C NPTIIC

**Immunogen:**

Recombinant Human Sodium-dependent phosphate transport protein 2C protein (1-76AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A