

# LANCL2 Antibody, Biotin conjugated



PACO58479

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Necessary for abscisic acid, (ABA) binding on the cell membrane and activation of the ABA signaling pathway in granulocytes.

**Gene ID:**

LANCL2

**Uniprot**

Q9NS86

**Synonyms:**

LanC-like protein 2 (Testis-specific adriamycin sensitivity protein), LANCL2, GPR69B TASP

**Immunogen:**

Recombinant Human LanC-like protein 2 protein (2-106AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A