

PTTG2 Antibody, FITC conjugated



PACO58542

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleus, peptidase inhibitor activity, homologous chromosome segregation, negative regulation of mitotic sister chromatid separation, regulation of transcription, DNA-templated.

Gene ID:

PTTG2

Uniprot

Q9NZH5

Synonyms:

Securin-2 (Pituitary tumor-transforming gene 2 protein), PTTG2

Immunogen:

Recombinant Human Securin-2 protein (1-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A