

PDP1 Antibody, FITC conjugated



PACO58546

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the dephosphorylation and concomitant reactivation of the alpha subunit of the E1 component of the pyruvate dehydrogenase complex.

Gene ID:

PDP1

Uniprot

Q9P0J1

Synonyms:

[Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 1, mitochondrial (PDP 1) (EC 3.1.3.43) (Protein phosphatase 2C) (Pyruvate dehydrogenase phosphatase catalytic subunit 1) (PDPC 1), PDP1, PDP PPM2C

Immunogen:

Recombinant Human [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 1, mitochondrial protein (1-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A