

FICD Antibody, Biotin conjugated



PACO58551

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Adenylyltransferase that mediates the addition of adenosine 5'-monophosphate (AMP) to specific residues of target proteins. Able to inactivate Rho GTPases in vitro by adding AMP to RhoA, Rac and Cdc42. It is however unclear whether it inactivates GTPases in vivo and physiological substrates probably remain to be identified.

Gene ID:

FICD

Uniprot

Q9BVA6

Synonyms:

Adenosine monophosphate-protein transferase FICD (EC 2.7.7. n1) (AMPylator FICD) (De-AMPylase FICD) (EC 3.1.4) (FIC domain-containing protein) (Huntingtin yeast partner E) (Huntingtin-interacting protein 13) (HIP-13) (Huntingtin-interacting protein E), FICD, HIP13 HYPE

Immunogen:

Recombinant Human Adenosine monophosphate-protein transferase FICD protein (63-185AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A