

DMRT3 Antibody, FITC conjugated



PACO58582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcription factor that plays a role in configuring the spinal circuits controlling stride in vertebrates. Involved in neuronal specification within specific subdivision of spinal cord neurons and in the development of a coordinated locomotor network controlling limb movements. May regulate transcription during sexual development.

Gene ID:

DMRT3

Uniprot

Q9NQL9

Synonyms:

Doublesex- and mab-3-related transcription factor 3, DMRT3, DMRTA3

Immunogen:

Recombinant Human Doublesex- and mab-3-related transcription factor 3 protein (92-224AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A