

# LARP6 Antibody, FITC conjugated



PACO58586

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Regulates the coordinated translation of type I collagen alpha-1 and alpha-2 mRNAs, CO1A1 and CO1A2. Stabilizes mRNAs through high-affinity binding of a stem-loop structure in their 5' UTR. This regulation requires VIM and MYH10 filaments, and the helicase DHX9.

**Gene ID:**

LARP6

**Uniprot**

Q9BRS8

**Synonyms:**

La-related protein 6 (Acheron) (Achn) (La ribonucleoprotein domain family member 6), LARP6

**Immunogen:**

Recombinant Human La-related protein 6 protein (294-453AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A