

# TIMELESS Antibody, FITC conjugated



PACO58606

---

## Product Information

|                               |  |
|-------------------------------|--|
| <b>Size:</b>                  | <b>Protein Background:</b>   |
| 50ug                          |  |
| <b>Reactivity:</b>            | <b>Gene ID:</b>  |
| Human                         | TIMELESS   |
| <b>Source:</b>                | <b>Uniprot</b>   |
| Rabbit                        | Q9UNS1   |
| <b>Isotype:</b>               | <b>Synonyms:</b>   |
| IgG                           | Protein timeless homolog (hTIM), TIMELESS, TIM TIM1 TIMELESS1                  |
| <b>Applications:</b>          | <b>Immunogen:</b>  |
| ELISA                         | Recombinant Human Protein timeless homolog protein (1004-1131AA).              |
| <b>Recommended dilutions:</b> | <b>Storage:</b>  |
|                               | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |

## Product Images

---

N/A

N/A