

TIMELESS Antibody, Biotin conjugated



PACO58607

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	TIMELESS
Source:	Uniprot
Rabbit	Q9UNS1
Isotype:	Synonyms:
IgG	Protein timeless homolog (hTIM), TIMELESS, TIM TIM1 TIMELESS1
Applications:	Immunogen:
ELISA	Recombinant Human Protein timeless homolog protein (1004-1131AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A