

EGR3 Antibody



PACO58624

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Probable transcription factor involved in muscle spindle development.

Gene ID:

EGR3

Uniprot

Q06889

Synonyms:

Early growth response protein 3 (EGR-3) (Zinc finger protein pilot), EGR3, PILOT


Immunogen:

Recombinant Human Early growth response protein 3 protein (201-281AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: MCF-7 whole cell lysate, A549 whole cell lysate. All lanes: EGR3 antibody at 3.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 43, 39 kDa. Observed band size: 43 kDa.