

ADAM23 Antibody, Biotin conjugated



PACO58680

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in cell-cell and cell-matrix interactions. This is a non-catalytic metalloprotease-like protein.

Gene ID:

ADAM23

Uniprot

O75077

Synonyms:

Disintegrin and metalloproteinase domain-containing protein 23 (ADAM 23) (Metalloproteinase-like, disintegrin-like, and cysteine-rich protein 3) (MDC-3), ADAM23, MDC3

Immunogen:

Recombinant Human Disintegrin and metalloproteinase domain-containing protein 23 protein (498-647AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A