

EMP2 Antibody, Biotin conjugated



PACO58740

Product Information

| | |
|-------------------------------|--|
| Size: | Protein Background: |
| 50ug | |
| Reactivity: | Gene ID: |
| Human | EMP2 |
| Source: | Uniprot |
| Rabbit | P54851 |
| Isotype: | Synonyms: |
| IgG | Epithelial membrane protein 2 (EMP-2) (Protein XMP), EMP2, XMP |
| Applications: | Immunogen: |
| ELISA | Recombinant Human Epithelial membrane protein 2 protein (22-66AA). |
| Recommended dilutions: | Storage: |
| | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |

Product Images

N/A

N/A