

NLRP6 Antibody, FITC conjugated



PACO58827

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NLRP6

Uniprot

P59044

Synonyms:

NACHT, LRR and PYD domains-containing protein 6 (Angiotensin II/vasopressin receptor) (PYRIN-containing APAF1-like protein 5), NLRP6, NALP6 PYPAF5

Immunogen:

Recombinant Human NACHT, LRR and PYD domains-containing protein 6 protein (154-306AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A