

PEX11A Antibody, HRP conjugated



PACO58834

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in peroxisomal proliferation and may regulate peroxisomes division. May mediate binding of coatomer proteins to the peroxisomal membrane. Promotes membrane protrusion and elongation on the peroxisomal surface.

Gene ID:

PEX11A

Uniprot

O75192

Synonyms:

Peroxisomal membrane protein 11A (HsPEX11p) (28 kDa peroxisomal integral membrane protein) (PMP28) (Peroxin-11A) (Peroxisomal biogenesis factor 11A) (Protein PEX11 homolog alpha) (PEX11-alpha), PEX11A, PEX11

Immunogen:

Recombinant Human Peroxisomal membrane protein 11A protein (106-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A