

# PEX11A Antibody, Biotin conjugated



PACO58836

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in peroxisomal proliferation and may regulate peroxisomes division. May mediate binding of coatomer proteins to the peroxisomal membrane. Promotes membrane protrusion and elongation on the peroxisomal surface.

**Gene ID:**

PEX11A

**Uniprot**

O75192

**Synonyms:**

Peroxisomal membrane protein 11A (HsPEX11p) (28 kDa peroxisomal integral membrane protein) (PMP28) (Peroxin-11A) (Peroxisomal biogenesis factor 11A) (Protein PEX11 homolog alpha) (PEX11-alpha), PEX11A, PEX11

**Immunogen:**

Recombinant Human Peroxisomal membrane protein 11A protein (106-219AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A