

PGLYRP3 Antibody, HRP conjugated



PACO58842

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Pattern receptor that binds to murein peptidoglycans (PGN) of Gram-positive bacteria. Has bactericidal activity towards Gram-positive bacteria. May kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. Binds also to Gram-negative bacteria, and has bacteriostatic activity towards Gram-negative bacteria. Plays a role in innate immunity.

Gene ID:

PGLYRP3

Uniprot

Q96LB9

Synonyms:

Peptidoglycan recognition protein 3 (Peptidoglycan recognition protein I-alpha) (PGLYRP1alpha) (PGRP-I-alpha) (Peptidoglycan recognition protein intermediate alpha), PGLYRP3, PGRPIA

Immunogen:

Recombinant Human Peptidoglycan recognition protein 3 protein (90-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A