

# PHACTR1 Antibody, FITC conjugated



PACO58847

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds actin monomers (G actin) and plays a role in the reorganization of the actin cytoskeleton and in formation of actin stress fibers. Plays a role in cell motility. Plays a role in the formation of tubules by endothelial cells. Regulates PPP1CA activity. Required for normal cell survival.

**Gene ID:**

PHACTR1

**Uniprot**

Q9C0D0

**Synonyms:**

Phosphatase and actin regulator 1, PHACTR1, KIAA1733 RPEL1

**Immunogen:**

Recombinant Human Phosphatase and actin regulator 1 protein (263-387AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A