

PREPL Antibody, Biotin conjugated



PACO58856

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable serine peptidase whose precise substrate specificity remains unclear. Does not cleave peptides after a arginine or lysine residue. May play a role in the regulation of synaptic vesicle exocytosis.

Gene ID:

PREPL

Uniprot

Q4J6C6

Synonyms:

Prolyl endopeptidase-like (EC 3.4.21) (Prolylendopeptidase-like), PREPL, KIAA0436

Immunogen:

Recombinant Human Prolyl endopeptidase-like protein (604-699AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A