

# RNASEH1 Antibody, FITC conjugated



PACO58875

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Endonuclease that specifically degrades the RNA of RNA-DNA hybrids. Plays a role in RNA polymerase II (RNAP II) transcription termination by degrading R-loop RNA-DNA hybrid formation at G-rich pause sites located downstream of the poly(A) site and behind the elongating RNAP II.

**Gene ID:**

RNASEH1

**Uniprot**

O60930

**Synonyms:**

Ribonuclease H1 (RNase H1) (EC 3.1.26.4) (Ribonuclease H type II), RNASEH1, RNH1

**Immunogen:**

Recombinant Human Ribonuclease H1 protein (73-185AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A