

RNASEH1 Antibody, Biotin conjugated



PACO58876

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Endonuclease that specifically degrades the RNA of RNA-DNA hybrids. Plays a role in RNA polymerase II (RNAP II) transcription termination by degrading R-loop RNA-DNA hybrid formation at G-rich pause sites located downstream of the poly(A) site and behind the elongating RNAP II.

Gene ID:

RNASEH1

Uniprot

O60930

Synonyms:

Ribonuclease H1 (RNase H1) (EC 3.1.26.4) (Ribonuclease H type II), RNASEH1, RNH1

Immunogen:

Recombinant Human Ribonuclease H1 protein (73-185AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A