

RPS23 Antibody, Biotin conjugated



PACO58880

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, cytosolic small ribosomal subunit, membrane, nucleoplasm, ribosome, RNA binding, structural constituent of ribosome, nuclear-transcribed mRNA catabolic process, nonsense-mediated decay, rRNA processing, SRP-dependent cotranslational protein targeting to membrane.

Gene ID:

RPS23

Uniprot

P62266

Synonyms:

40S ribosomal protein S23 (Small ribosomal subunit protein uS12), RPS23

Immunogen:

Recombinant Human 40S ribosomal protein S23 protein (2-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A