

SASH1 Antibody



PACO58889

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:200-1:500

Protein Background:

May have a role in a signaling pathway. Could act as a tumor suppressor.

Gene ID:

SASH1

Uniprot

O94885

Synonyms:

SAM and SH3 domain-containing protein 1 (Proline-glutamate repeat-containing protein), SASH1, KIAA0790 PEPE1


Immunogen:

Recombinant Human SAM and SH3 domain-containing protein 1 protein (76-191AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



IHC image of PACO58889 diluted at 1:300 and staining in paraffin-embedded human placenta tissue performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.