

# TPK1 Antibody, HRP conjugated



PACO58954

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the phosphorylation of thiamine to thiamine pyrophosphate. Can also catalyze the phosphorylation of pyrithiamine to pyrithiamine pyrophosphate.

**Gene ID:**

TPK1

**Uniprot**

Q9H3S4

**Synonyms:**

Thiamin pyrophosphokinase 1 (hTPK1) (EC 2.7.6.2) (Placental protein 20) (PP20) (Thiamine pyrophosphokinase 1), TPK1

**Immunogen:**

Recombinant Human Thiamin pyrophosphokinase 1 protein (1-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A