

ZNF3 Antibody, FITC conjugated



PACO58979

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in cell differentiation and/or proliferation.

Gene ID:

ZNF3

Uniprot

P17036

Synonyms:

Zinc finger protein 3 (Zinc finger protein HF.12) (Zinc finger protein HZF3.1) (Zinc finger protein KOX25), ZNF3, KOX25

Immunogen:

Recombinant Human Zinc finger protein 3 protein (92-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A