

# USB1 Antibody, HRP conjugated



PACO59010

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Phosphodiesterase responsible for the U6 snRNA 3' end processing. Acts as an exoribonuclease (RNase) responsible for trimming the poly(U) tract of the last nucleotides in the pre-U6 snRNA molecule, leading to the formation of mature U6 snRNA 3' end-terminated with a 2',3'-cyclic phosphate.

**Gene ID:**

USB1

**Uniprot**

Q9BQ65

**Synonyms:**

U6 snRNA phosphodiesterase (hUsb1) (EC 3.1.4), USB1, C16orf57

**Immunogen:**

Recombinant Human U6 snRNA phosphodiesterase protein (1-90AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A