

HES5 Antibody, Biotin conjugated



PACO59132

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional repressor of genes that require a bHLH protein for their transcription. Plays an important role as neurogenesis negative regulator.

Gene ID:

HES5

Uniprot

Q5TA89

Synonyms:

Transcription factor HES-5 (Class B basic helix-loop-helix protein 38) (bHLHb38) (Hairy and enhancer of split 5), HES5, BHLHB38

Immunogen:

Recombinant Human Transcription factor HES-5 protein (1-166AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A