

SIRPB1 Antibody, HRP conjugated



PACO59134

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Immunoglobulin-like cell surface receptor involved in the negative regulation of receptor tyrosine kinase-coupled signaling processes.

Gene ID:

SIRPB1

Uniprot

Q5TFQ8

Synonyms:

Signal-regulatory protein beta-1 isoform 3 (SIRP-beta-1 isoform 3), SIRPB1

Immunogen:

Recombinant Human Signal-regulatory protein beta-1 isoform 3 protein (30-206AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A