

# ULK2 Antibody, FITC conjugated



PACO59223

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ULK2

**Uniprot**

Q8IYT8

**Synonyms:**

Serine/threonine-protein kinase ULK2 (EC 2.7.11.1) (Unc-51-like kinase 2), ULK2, KIAA0623

**Immunogen:**

Recombinant Human Serine/threonine-protein kinase ULK2 protein (211-276AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A