

GPT2 Antibody, Biotin conjugated



PACO59248

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the reversible transamination between alanine and 2-oxoglutarate to form pyruvate and glutamate.

Gene ID:

GPT2

Uniprot

Q8TD30

Synonyms:

Alanine aminotransferase 2 (ALT2) (EC 2.6.1.2) (Glutamate pyruvate transaminase 2) (GPT 2) (Glutamic--alanine transaminase 2) (Glutamic--pyruvic transaminase 2), GPT2, AAT2 ALT2

Immunogen:

Recombinant Human Alanine aminotransferase 2 protein (95-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A