

SEMA6D Antibody, Biotin conjugated



PACO59316

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Shows growth cone collapsing activity on dorsal root ganglion (DRG) neurons in vitro. May be a stop signal for the DRG neurons in their target areas, and possibly also for other neurons. May also be involved in the maintenance and remodeling of neuronal connections.

Gene ID:

SEMA6D

Uniprot

Q8NFY4

Synonyms:

Semaphorin-6D, SEMA6D, KIAA1479

Immunogen:

Recombinant Human Semaphorin-6D protein (913-1068AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A