

MRAP Antibody, HRP conjugated



PACO59334

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Modulator of melanocortin receptors (MC1R, MC2R, MC3R, MC4R and MC5R). Acts by increasing ligand-sensitivity of melanocortin receptors and enhancing generation of cAMP by the receptors. Required both for MC2R trafficking to the cell surface of adrenal cells and for signaling in response to corticotropin (ACTH). May be involved in the intracellular trafficking pathways in adipocyte cells.

Gene ID:

MRAP

Uniprot

Q8TCY5

Synonyms:

Melanocortin-2 receptor accessory protein (B27) (Fat cell-specific low molecular weight protein) (Fat tissue-specific low MW protein), MRAP, C21orf61 FALP

Immunogen:

Recombinant Human Melanocortin-2 receptor accessory protein (59-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A