

# SEMA3C Antibody, FITC conjugated



PACO59367

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds to plexin family members and plays an important role in the regulation of developmental processes. Required for normal cardiovascular development during embryogenesis. Functions as attractant for growing axons, and thereby plays an important role in axon growth and axon guidance.

**Gene ID:**

SEMA3C

**Uniprot**

Q99985

**Synonyms:**

Semaphorin-3C (Semaphorin-E) (Sema E), SEMA3C, SEMAE

**Immunogen:**

Recombinant Human Semaphorin-3C protein (292-396AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A