

# SLC19A3 Antibody, HRP conjugated



PACO59378

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Mediates high affinity thiamine uptake, probably via a proton anti-port mechanism. Has no folate transport activity.

**Gene ID:**

SLC19A3

**Uniprot**

Q9BZV2

**Synonyms:**

Thiamine transporter 2 (ThTr-2) (ThTr2) (Solute carrier family 19 member 3), SLC19A3

**Immunogen:**

Recombinant Human Thiamine transporter 2 protein (191-282AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A