

SLC19A3 Antibody, Biotin conjugated



PACO59380

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates high affinity thiamine uptake, probably via a proton anti-port mechanism. Has no folate transport activity.

Gene ID:

SLC19A3

Uniprot

Q9BZV2

Synonyms:

Thiamine transporter 2 (ThTr-2) (ThTr2) (Solute carrier family 19 member 3), SLC19A3

Immunogen:

Recombinant Human Thiamine transporter 2 protein (191-282AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A