

CNNM2 Antibody, HRP conjugated



PACO59410

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Divalent metal cation transporter. Mediates transport of divalent metal cations in an order of Mg(2+) > Co(2+) > Mn(2+) > Sr(2+) > Ba(2+) > Cu(2+) > Fe(2+).

Gene ID:

CNNM2

Uniprot

Q9H8M5

Synonyms:

Metal transporter CNNM2 (Ancient conserved domain-containing protein 2) (Cyclin-M2), CNNM2, ACDP2

Immunogen:

Recombinant Human Metal transporter CNNM2 protein (739-875AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A