

ALX4 Antibody, FITC conjugated



PACO59487

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcription factor involved in skull and limb development. Plays an essential role in craniofacial development, skin and hair follicle development.

Gene ID:

ALX4

Uniprot

Q9H161

Synonyms:

Homeobox protein aristaless-like 4, ALX4, KIAA1788

Immunogen:

Recombinant Human Homeobox protein aristaless-like 4 protein (112-216AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A